

# REFURBISHED MOLLER EOS 900 ON 3011 FLOORSTAND

## MOLLER OPTICS

- Sharp Optical View
- Terrific Depth of Focus
- Widefield 0-75 Degree Binocular
- High Resolution, Internal, Motorized Focus and Zoom with Re-centering
- Three Angles of Illumination: Vertical for Bright Red Reflex, Coaxial for Contrast and Oblique for Extraocular Procedures
- 1:6 Magnification Ratio
- Automatic Speed Control for X/Y and Focus
- Fully Illuminated Visual Field with Manual Spot Light

## MOLLER 3011 FLOORSTAND

- Stable, Movable Base
- Long, Articulating Arm with Smart Brakes in Each Joint and Internal Duct for All Cables
- Smart Electronics with Self-Check of All Control Functions
- Independent Digital Speed Controls for Focus, Zoom & XY
- Adjustable Spring Tension for Balancing
- Durable Friction Knobs for Rotation of Arms and Up/Down Articulation
- Locking Spindles to Secure or Level Base
- Easy Maneuvering with Side Handles
- Sturdy, Oversized Casters

## FIBEROPTIC ILLUMINATOR

- Two Bulb System for Immediate Back-up
- Varying Rheostat for High to Low Illumination Levels
- Reliable and Easy to Exchange Bulbs
- Selector for 110V or 220V

## XY FUNCTION

- Planar or Rotational Movement
- Automatic Centering Button
- Motorized Movement Activated by Joystick on Foot Control

## FOOT CONTROL

- Push Button for Illumination On/Off Available
- Smooth Pedals for Motorized Focus and Zoom
- Movable Joystick for Motorized XY
- Durable, Water-Resistant Rubber



Complete  
2 Year  
Warranty

# ENDURE MEDICAL, INC.

## Optional Accessories



The 2011 floorstand, various assistant scopes, footpedals and replacement parts are always available.



New Video Systems and Accessories for All Microscopes



We are a family business that specializes in new and refurbished surgical microscopes and accessories. We've been in business since 1986 and offer a complete 2-year warranty (including parts and labor) on every system. Our goals include providing superior products, competitive pricing and consistent customer service. We promise to conduct business in a timely, accurate and professional manner. Please call on us to discover what a small business of quality individuals can do for you!

## Technical Specifications

Microscope	Specifications
Binocular	0-75 Inclinable ; f=170mm
Objective Lens	175mm, 200mm or your choice
Magnification	Motorized 1:6, .4...2. 4
Focus	Motorized with $\pm 30$ mm range
Eyepieces	10X 22B Widefield
Pupillary Distance	55mm - 75mm
Filter	Standard UV
Illumination	
Lightsource	Integrated 12V 75W
Illuminated Field	56 - 32mm
Power Source	110 - 240V ; 50 - 60Hz
Mounting System	
Type	FS 3011 ; FS 2011 Also Available
Horizontal Reach	Arm = 52" ; End of Optics = 54"
Vertical Reach	$\pm 72$ "
Rotation Angle	Lamp Arm = 360 degrees Articulated Arm = 270 degrees
Weight	
Shipping Weight	$\pm 440$ lbs
Warranty	
Term	2 years including parts and labor

1455 Ventura Drive  
Cumming, GA 30040  
[www.enduremed.com](http://www.enduremed.com)



Phone: 800-736-3873  
Fax: 770-888-3991  
Local: 770-888-3755